

## North2South Chemiluminescent Detection

Fisher catalog # 17097

<https://www.thermofisher.com/order/catalog/product/17097>

Danna Eickbush

Note: *Volumes may change depending on size of membrane & containers used.*  
Rinse forceps between each step.

1. The day before (ie after starting hybridization) bring Hybridization Stringency Wash Buffer (2X) to room temperature- warm to 37°C if there is still a precipitate. Then make dilution for needed volume. For 1 blot: 40 mls of 2X buffer plus 40 mls of sterile H<sub>2</sub>O. This wash buffer can be left in a water bath overnight at the hybridization temperature. Note: *Stringency of wash buffer can be adjusted (0.1X to 2X).*
2. The next day, put Blocking Buffer, 4X Blocking Wash Buffer, & sterile water in a 37°C water bath. Put Substrate Equilibration Buffer, Peroxide Solution, & Luminol/Enhancer Solution at room temperature.
3. Pour off hybridization solution. Wash membrane with "1X" Hyb washing buffer.
  - a. 25 mls, 15-20 minutes, hyb temp, shaking or rotating
  - b. 25 mls, 15-20 minutes, hyb temp, shaking or rotating
  - c. 25 mls, 15-20 minutes, hyb temp, shaking or rotating
4. Put membrane in a clean box (slightly larger than the membrane).
  - a. Add 16 mls of Blocking Buffer.
  - b. Shake gently, room temperature, 15 minutes
5. During step 4, mix 16 mls of Blocking Buffer + 55 ul of the strep/HRP conjugate.
  - a. After the 15 minute incubation in step 4, decant blocking buffer.
  - b. Add Blocking Buffer + conjugate.
  - c. Shake gently, room temperature, 15 minutes.
6. During step 5, make 1X Blocking Wash Buffer. For 1 membrane, mix 25 mls of prewarmed 4X Wash buffer and 75 mls of pre-warmed H<sub>2</sub>O. Keep 1X buffer at 37°C.
  - a. After the 15 minute incubation in step 5, transfer membrane to a clean box.
  - b. Rinse briefly with 20 mls of 1X block wash.
  - c. Wash membrane 20 mls, 5 minutes, gentle shaking, room temp
  - d. Wash membrane 20 mls, 5 minutes, gentle shaking, room temp
  - e. Wash membrane 20 mls, 5 minutes, gentle shaking, room temp
  - f. Wash membrane 20 mls, 5 minutes, gentle shaking, room temp

7. Transfer membrane to a clean box.
  - a. Wash membrane 5 minutes in Substrate Equilibration Buffer, gentle shaking, room temperature
  
8. Mix 6 mls luminol/enhancer and 6 mls peroxide solutions. Protect from light during this and following steps.
  - a. Pour mixture into clean container
  - b. Using clean forceps, pick up membrane, blot one corner on filter paper for 3 seconds, place membrane DNA/RNA side down onto the luminol/proxide "puddle".
  - c. Incubate 5 minutes, no shaking, room temperature.
  - d. Pick up membrane with forceps, blot one corner on filter paper for 3 seconds, wrap in saran wrap.
  
9. Expose membrane to CCD camera. Adjust exposure length as needed.