

RNA Isolation

Danna Eickbush

Special Cautions:

Steps 3 and 6

Phenol (corrosive, toxic), chloroform/sevag (anesthetic, toxic)- work in hood, wear proper PPE, dispose waste in proper hazardous waste containers.

Step 8

Ethidium bromide (mutagen)- wear gloves, dispose waste in proper hazardous waste containers. UV radiation (skin burns, eye damage)- cover all exposed skin, wear eye goggles.

1. 25 females or 40 males.
2. Grind flies in 200 ul Buffer A.
Add 200 ul Buffer B and 10 ul proteinase K (10 mg/ml).

Incubate 1.5 hours at 37⁰C. Add another 10 ul of proteinase K and continue incubation for another 1.5 hours.
3. Phenol/sevag extract (400 ul) 2X.
Sevag extract (400 ul) 1X.

Add 40 ul 3M NaAcetate and 900 ul cold 95% ethanol.
Gently mix. Place at -20⁰C for 25 minutes. Centrifuge 30 minutes, 13K in coldroom.
4. Dump off ethanol, rinse with 70% ethanol, quick spin, pipette off excess ethanol. Air dry for a few minutes. DO NOT OVER DRY.
5. Resuspend in 200 ul DNase I buffer (40 mM Tris 7.9, 5 mM MgCl₂, 5 mM CaCl₂).
Add 20 units of DNase I (Ambion). Incubate at 37⁰C for 40 minutes.
6. Phenol/sevag extract (200 ul) 2X.
Sevag extract (200 ul) 1X.

Add 20 ul 3M NaAcetate, 600 ul ethanol. -20⁰C overnight. Centrifuge 25 minutes in cold. Rinse pellet with 70% ethanol, quick spin, pipette off excess ethanol, air dry a few minutes.
7. Resuspend 25 ul water.
8. Determine concentration, check integrity on agarose/ethidium bromide gel.
The yield is usually between 75 and 100 ug total RNA.

Buffer A

10 mM Tris 8.0
60 mM NaCl
10 mM EDTA
0.15 mM spermidine
0.15 mM spermine
5% sucrose

Buffer B

200 mM Tris 9.0
30 mM EDTA
2% SDS
5% sucrose